



VUNO Med®-Fundus AI™

VUNO Med® Fundus AI™

AI-based Screening Solution
for Various Abnormalities in the Fundus



Seoul, South Korea (HQ) 9F, 479, Gangnam-daero, Seocho-gu, Seoul, Korea
T +82.2.515.6646 | F +82.2.515.6647 | www.vuno.co





VUNO Med®-Fundus AI™

VIEW THE INVISIBLE
KNOW THE UNKNOWN



Impacts of Vision Impairment

Globally, the leading causes of vision impairment and blindness include age-related macular degeneration (AMD), diabetic retinopathy (DR), glaucoma and cataract.

170 million individuals are affected by AMD globally, which affects the central vision and impacts quality of life making carrying out daily activities difficult

415 millions adults globally were living with diabetes in 2015 and it is estimated that approximately 93 million people worldwide have DR and one third of diabetic population will develop some form of retinopathy in their lifetime

It is estimated that 57.5 million people are affected by primary open-angle glaucoma (POAG) worldwide and by 2020 the numbers will rise to 76 million and further rising up to 111.8 million by 2040

Governments around the world are bannng together to address the global magnitude of preventable blindness and implement policies to address and manage visual impairments. Early intervention, effective and accessible eye care services, like eye screening, is key for effectively preventing visual impairments and vision loss.

Eye screening is being integrated as a core part of overall eye and/or disease management and care. Fundus examinations allow early detection of pre-symptomatic eye disease and subsequently allows patients to receive early treatments before vision loss becomes irreversible.

Effective Fundus Screening is Not Easy

Shortage of ophthalmologists

There is not enough ophthalmologists to treat and perform fundus screening for increasing number of patients due to population growth and aging

Challenge for ophthalmologists to reach high risk patients

Diabetic and hypertensive patients who are at risk of developing eye complications visit primary care clinics or GP first before getting referred to ophthalmologists

Lack of screening tool and knowledge in primary care settings

Primary care practices usually are not equipped with eye screening tool (in this case, fundus camera) and image interpretation knowledge


Unnecessary referrals lead to increased workload and cost

Without proper knowledge and consultation in primary care settings, even normal patients are sometimes referred to ophthalmology

Interrupted and ineffective workflow in telemedicine service

Telemedicine or teleophthalmology service needs to send fundus image for interpretation by trained grader before they can present it to patients

View More with
VUNO Med[®]
Fundus AI[™]




What abnormalities do you see?



- | | |
|---|---|
| <input type="checkbox"/> Hemorrhage | <input type="checkbox"/> Non-glaucomatous Disc Change |
| <input type="checkbox"/> Hard Exudate | <input type="checkbox"/> RNFL Defect |
| <input type="checkbox"/> Cotton Wool Patch | <input checked="" type="checkbox"/> Membrane |
| <input checked="" type="checkbox"/> Drusen & Drusenoid Deposits | <input type="checkbox"/> Chorioretinal Atrophy / Scar |
| <input type="checkbox"/> Vascular Abnormality | <input checked="" type="checkbox"/> Macular Hole |
| <input type="checkbox"/> Glaucomatous Disc Change | <input type="checkbox"/> Myelinated Nerve Fiber |

What does VUNO Med®-Fundus AI™ do?

VUNO Med®-Fundus AI™ is an AI-based screening solution that provides diagnostic support for common ocular diseases based on the detection of 12 retinal findings in the fundus of the eye.



12

Detection of 12 retinal findings associated with diagnosis of eye diseases



95%

High sensitivity and specificity with AUROC 96.2 – 99.9%



2s

Immediate fundus interpretation in only 2 seconds



100,000+

Developed and validated with 100,000+ images graded by 57 ophthalmologists of different subspecialties

Vision-Threatening Eye Diseases



Diabetic Retinopathy



Age-related Macular Degeneration



Retinal Vein Occlusion



Glaucoma Suspect



Epiretinal Membrane



Macular Hole

How can I use
VUNO Med[®]-Fundus AI[™]
to optimize and improve
my clinical workflow?



Diagnostic Accuracy

- Accurately detects and locates (AUROC = 95% and above) 12 fundus abnormalities trained on 103,262 fundus images with 57 Ophthalmologists
- Provides as a cross check reference and increase confidence in interpretation among untrained graders
- Improves sensitivity – ensure no findings are missed out

Effective Screening

- Provides interpretation results within 2 seconds
- Provides more effective and necessary referrals to high volume ophthalmology clinics, so patient consultation and treatment can be prioritized
- Provides tangible evidence to patients for referral

Consistency

- Provides consistent diagnosis results without depending on interpreter's knowledge and experience



VUNO Med®-Fundus AI™ is Intuitive & User-friendly



1 Step

Take and Upload Fundus Image

A high-quality fundus image taken with mydriatic, or non-mydriatic fundus camera is uploaded in **JPG, DCM or PNG format**



2 Step

Perform Analysis with VUNO Med®-Fundus AI™

In only **2 seconds**, VUNO Med®-Fundus AI™ algorithm provides **automatic detection** of abnormal findings




3 Step

Receive Examination Report

Report is generated with labeled abnormal findings for the examiner and/or patients for further consultation

VUNO Med®-Fundus AI™ Report

Patient report containing visualized abnormality findings that facilitate doctor-patient communication and enhance patient satisfaction.



VUNO Med®-Fundus AI™ is Versatile & Easy-to-Implement



Cloud-Based

- A virtual server where **VUNO Med®-Fundus AI™** is installed and customers can access the solution by logging on to the server (URL) to use our products
- **VUNO Med®-Fundus AI™** can be used on a license / subscription basis or pay-per-scan basis



On-Premise

- A virtual server where **VUNO Med®-Fundus AI™** is installed and run within a client's in-house infrastructure
- **VUNO Med®-Fundus AI™** can be used on a license / subscription basis or pay-per-scan basis



On-Device

- A virtual server where **VUNO Med®-Fundus AI™** is installed and embedded in diagnostic modalities through the SDK integration
- **VUNO Med®-Fundus AI™** can be used as a bundled product in diagnostic modality



How and where can VUNO Med[®]-Fundus AI[™] be used?

Ophthalmology

- Diagnostic support tool for eye hospitals and clinics' eye screening campaigns or programs

Primary Care Settings & Screening Centers

- Diagnostic support tool for common eye diseases for effective patient referral to ophthalmology
- Offer as additional ophthalmology screening service for patients

Tele-Ophthalmology, Tele-Medicine and Medical Imaging Centers

- Time and cost saving diagnostic support tool
- Offer as additional ophthalmology screening service

NGOs & IGOs

- Effective screening tool and partnerships for eye screening campaigns or care programs

Camera Suppliers & Manufacturers

- Offer as competitive advantage over competitors
- New market penetration into non-ophthalmology settings

VUNO Med[®]-Fundus AI[™] is Versatile & Easy-to-Implement

1 Choose type of service



2 Choose payment type



3 Perform analysis at your fingertips



Clinical Data of VUNO Med[®]-Fundus AI[™]

VUNO Med®-Fundus AI™ is Comprehensively Trained & Validated


Developed and validated with 100,000+ images graded
by 57 ophthalmologists of different subspecialties

- Consistent and stable performance on internal test set (AUROC 96.2 – 99.9%)
- Consistent and stable performance on external test set (France and India) (AUROC 94.7 – 98.0%)

Diagnostic Performance

■ Hospital A (In-House test, Korea)
 ■ Hospital B (External test, India)
 ■ Hospital C (External test, France)

Performance Comparison (Hemorrhage case)



Publications

No	Year	Jour/Cof	Type	Title
1	2018	JKMS	Journal	A Novel Fundus Image Reading Tool for Efficient Generation of a Multi-dimensional categorical Image Database for Machine Learning Algorithm Training
2	2018	MICCAI	Workshop	Classification of Findings with Localized Lesions in Fundoscopic Images using a Regionally Guided CNN
3	2018	MICCAI	Workshop	An Efficient and Comprehensive Labeling Tool for Large-Scale Annotation of Fundus Images
4	2018	JDI	Journal	Towards Accurate Segmentation of Retinal Vessels and the Optic Disc in Fundoscopic Images with Generative Adversarial Networks
5	2018	JDI	Journal	Laterality Classification of Fundus Images Using Interpretable Deep Neural Network
6	2019	MIA	Journal	REFUGE Challenge: A unified framework for evaluating automated methods for glaucoma assessment from fundus photographs
7	2019	MIA	Journal	IDRiD: Diabetic Retinopathy – Segmentation and Grading Challenge
8	2019	Ophthalmology	Journal	Development and Validation of Deep Learning Models for Screening Multiple Abnormal Findings in Retinal Fundus Images

Awards



Korea's 1st Innovative Medical Device



1st place in Automatic Detection Challenge on Age-related Macular Degeneration, 1st place in Diabetic Retinopathy Segmentation and Grading Challenge



1st place in Retinal Fundus Glaucoma Challenge



IDRiD Challenge
All sub-challenges 1st and 2nd place



REFUGE Challenge
Glaucoma classification 1st place
Overall 2nd place



ODIR Challenge
Ocular Disease Intelligent Recognition
Top 2% (offsite)



PALM challenge
Pathological Myopia classification
1st place (off-site)



AGE Challenge
Angle closure classification
1st place (off-site)



DeepDRiD Challenge
DR grading (macular, disc-centered):
3rd place
DR grading (ultra-wide): 1st place



ADAM Challenge
Overall 1st place All sub-challenges 1st, 2nd places
except for disc segmentation (5th place)

VUNO Med®-Fundus AI™ Demo



Are you interested in trying out VUNO Med®-Fundus AI™?

Do you

Want to request a demo?

Want to learn more about how this technology improves and optimizes clinical workflows?

Want to discover how you can implement this technology?

Just let us know!





VIEW THE INVISIBLE KNOW THE UNKNOWN

VUNO is one of the biggest healthcare AI companies in the world with strong commercial traction of road-tested and clinically-proven solutions.

Founded back in December 2014, the company successfully obtained the 1st Korean FDA (MFDS) approval for medical AI software. Since then, VUNO has rapidly expanded its product portfolio to a diverse and comprehensive line up, that is reshaping the delivery of medical imaging diagnostics in the field of radiology, pathology, bio-signal and medical speech recognition.

A proven track record of success in winning international technical challenges coupled with a series of regulatory approval and commercial deployments to over 500 hospitals, make VUNO an undisputed leader in the global medical AI arena.



If you are interested in distributing **VUNO Med®-Fundus AI™** contact us.

Office

Seoul, South Korea (HQ)

9F, 479, Gangnam-daero, Seocho-gu, Seoul, Korea

T +82.2.515.6646 **F** +82.2.515.6647

Boston, USA (Branch Office)

Contact Us

Please feel free to email us about any inquiries or questions

sales@vuno.co

www.vuno.co

